

# Medical Suction Filters

FOR COMPLIANCE WITH SAFETY  
STANDARDS



## Comprehensive Range

- Odor filters, bacterial filters
- Virus filters new in 2018
- Medela and conical connections
- Disposable and reusable

## High-Performance

- Dedicated, high-performance medical membranes
- Lab-tested high bacteria retention grades
- Low air resistance due to the large surface area



ISO 10079-1/-3

**Medela offers a vast portfolio of disposable filters that are seamlessly compatible with all Medela healthcare devices. They enable you to custom-configure your suction system to suit a variety of hygienic requirements.**

- Bacteria filters provide the first, fundamental level of hygienic protection. A basic requirement if your system does not include an embedded bacterial filter.
- Overflow filters offer an essential safety mechanism to prevent fluids from entering the pump or central vacuum system. Two redundant overflow mechanisms should be present in any aspiration system.
- Odor membranes effectively trap gases and odors, especially useful in containers used over multiple hours.
- In addition, Medela offers reusable bacteria filters for applications in which membranes need to be replaced with higher frequency.
- To ensure safety for patients and hospital staff, the pump filter exchange frequency should be performed according to the instructions for use and in-line with applicable local infection control policy.



Medela Connection



Conical Connection

	Medela/ Medela Connections	Conic/ Conic Connections	Conic/ Medela Connections	Recommended application	For following collection systems
<b>Bacteria</b>	077.0573	077.0574		Basic hygienic protection	Reusable jars
<b>Bacteria + Overflow</b>	077.0571	077.0572	077.0834 (with tubing) 077.0577	Second stage overflow protection	Disposable liners Reusable jars
<b>Odor + Bacteria</b>	077.0575	077.0576		Reduction of odor	Disposable liners Reusable jars
<b>Reusable bacteria filter</b>	077.0581 (filter housing) 077.1014 (filter paper)			Alternative to disposable bacteria filters, especially for appl. with frequent filter exchange	Reusable jars

#### Filter exchange

- In case of overflow or clogging, or
- Weekly with single-patient use, or
- Daily with multiple patient use.

To detect a clogged filter, perform a filter test before each patient use. Comply with in-house directives on hygiene.



Medical Vacuum Technology for  
Healthcare Professionals

**Please contact us or your local  
Medela Representative for details.**

 Medela AG  
Lättichstrasse 4b  
6341 Baar, Switzerland  
www.medela.com

Local contact